

TeleFipho^o


Wireless Fiberphotometry

Fiber photometry is a powerful technique to detect calcium signal from specific neuron in awake animals using calcium indicator protein represented by GCaMP. In freely moving condition, the long optic fiber attached to the head of the animal can interfere with experimental setup, therefore can be a limitation factor of your project.

Our innovative new product, TeleFipho, includes all required components for fiber photometry - optic fiber, filter cube, light source, photo detector - and also wireless transmission circuit, in the very small 3g body. TeleFipho definitely will not block free behavior of your animals, enabling novel experimental approach using fiber photometry.

○ Features

- ⊗ World first commercial wireless fiber photometry
- ⊗ Small headstage / good for mice, rats, marmosets, etc.
- ⊗ Standard 2.5mm ferrule cannula
- ⊗ Rechargeable by a dedicated charger
- ⊗ Adjustable excitation LED power
- ⊗ Adjustable signal offset
- ⊗ For GCaMP or GFP-like indicators
e.g. dLight, GRAB-DA, GRAB-5HT, etc.



Specification	
Headstage weight	3g
Headstage size	12 x 12 x 22 mm
Excitation wavelength	LED peak 470nm, Filter band 445~490nm
Emission wavelength	Filter band 500~550nm
Excitation power	10~300μW@Fiber end (Adjustable)
Sampling rate	100Hz
AD resolution	16bit
Photo sensor	Photo diode
Gain	10 ¹⁰ V/A
Battery life	2hours@Excitation power 30μW
Transmission band	2.4GHz
Transmission distance	2m
Power	Battery powered, rechargeable
Receiver I/O	1x Photometry analog out, 1x General purpose analog In (0~5V)
PC Interface	USB / TeleFipho software (for Windows10)
Cannula	core: 400μm/NA0.39, Cladding: 425μm, Ferrule: 2.5mm core: 200μm/NA0.39, Cladding: 225μm, Ferrule: 2.5mm core: 600μm/NA0.39, Cladding: 630μm, Ferrule: 2.5mm


○ TeleFipho Standard Set

TeleFipho-set consists of the following items:


- TeleFiT** 1x TeleFipho Transmitter Headstage
- TeleFiR** 1x TeleFipho Receiver
- TeleFiCharger** 1x TeleFipho Charger
- TeleFiC-I-d** (※) 3x TeleFipho Cannula
- TeleFiTool** 1x Insertion Tool for TeleFipho
- TeleFiDummy** 1x Dummy Headstage
1x TeleFipho Software Installer

(※) Please specify cannula length.


5mm length and 400μm cannula will come without specifying.
e.g. Length 3.5mm, dia. 400μm: TeleFiC-3.5-400 (in mm, resolution 0.1mm)



○ Data



Activity synchronized GCaMP signal in mouse dorsal striatum.
By courtesy of Dr. Keitaro Yoshida in Kenji Tanaka Lab, Keio University.




Stress induced GCaMP signal from orexin neuron in mouse hypothalamus.
By courtesy of Dr. Daisuke Ono in Akihiro Yamanaka Lab, Nagoya University.

○ Stand Alone / External Recorder

TeleFipho receiver can be directly connected to PC via USB for fiber photometry data recording in TeleFipho software. TeleFipho software is good enough for simple raw data / low pass filtered data recording. Data can be directly read from Matlab / Octave, or can be exported to ASCII format directly readable by pMAT, open source analysis software for fiber photometry.

For more complicated online analysis and / or parallel recording of other physiological data, you can access photometry data as analog signal via BNC on the front panel of TeleFipho receiver so that you can connect external data recording system (e.g. PowerLab from ADI).



Good for mouse!

Model	Description
TeleFipho-set	TeleFipho Standard Set
TeleFiT	TeleFipho Transmitter Headstage
TeleFiR	TeleFipho Receiver
TeleFiCharger	TeleFipho Charger
TeleFiC-I-d	TeleFipho Cannula
TeleFiTool	Insertion Tool for TeleFipho
TeleFiDummy	TeleFipho Dummy Headstage

In US & Canada:

AMUZA INC

10060 Caroll Canyon Road, Suite 100, San Diego, California, USA, 92131
URL: <https://amuzainc.com>
Tel: (858) 225-6869 Fax: (858) 560-8040

Other Countries:



Bio Research Center

Towa-Takaoka Bldg. 4F, 2-28-24 Izumi, Higashi-ku, Nagoya, Japan 461-0001
URL: <http://www.teleopto.com> Mail: sales-intl@brck.co.jp
Tel: +81-52-932-6421 Fax: +81-52-932-6755